# **NGT** CLASSIC

## NORMALAB GUM TESTER



## **STANDARDS**

ASTM D 381, IP 131, IP 540, ISO 6246, DIN 51784 and related methods.

## SCOPE

This test method covers the determination of the existent gum content of aviation fuels, and the gum content of motor gasolines or other volatile distillates in their finished form at the time of test.



## **SPECIFICATIONS**

The device determines existent gum in fuels by steam-jet and air-jet evaporation.

### **Key features**

- Temperature range up to 260°C
- Switch air / steam
- Small footprint
- 5 positions
- Electronic temperature control
- Easy installation of the 5 injectors
- Integrated superheater for steam method

#### Safety

Powercut in case of overheat

## Quality

- Calibration kit for air
- Calibration kit for vapour

#### Accessories

- REF 9413080 Air compressor: 2 stages suitable to deliver correct flow of air, 300 I capacity. For use on AC 400 V - 3 ph - 7 A - (1670x600x760mm) - 125 kg
- REF 941309 **Steam generator**: (400V 3 phases) for use with ref. 941320 apparatus, equipped with sensitive device to maintain vapour output at requested quantity, automatic water feeding by means of a water pump, automatic water level control, two safety devices (pressure switch and safety valve), manometer, liquid level and supplied with flexible connecting pipe. *Supply*: 400V, III, 18A, 50-60 Hz 12 kW (600x450x1100mm) 75 kg
- REF 41520 Analytical balance: 120 grams (0,1 mg)

## ORDERING INFORMATION

941320

## **Scope of delivery:**

NGT Classic is delivered ready to use with:

- 1 ASTM Thermometer 3C (P/N 11491)\*
- 5 Beakers (P/N 19035)
- 5 Conical outlets (P/N 941306)
- 1 Handling equipment (Forceps) (P/N 10686)
- 1 Pair of gloves (P/N 9417444)

\*Without certificate, calibration on request

#### Site requirements:

- Power supply: AC 230 V, 50/60 Hz 16 A
- Dimension:(W) 580x (D) 350x (H) 440 mm
  (W) 580x (D) 350x (H) 630 mm with gripp
- Weight: 40 kg

Required pressure: 2,15 Bars max

According to your application, accessories and options must be ordered separately.



**CONTACT**: sales@normalab.com

Normalab FRANCE SAS ZA Caux Multipôles 1 - 76190 Valliquerville

Tel.: +33 232.700.100 Fax: +33 232.704.732 **DISTRIBUTED BY**